



ABHD10 Antibody

Product Code	CSB-PA591994
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NUJ1
Immunogen	Synthesized peptide derived from human protein . at AA range: 220-300
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	WB ,ELISA; WB:1:500-2000, ELISA:1:5000-20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ABHD10