



# S Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA33245YB11GMY
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P0DTC2
<b>Immunogen</b>	Recombinant Human Novel Coronavirus Spike glycoprotein (S) (16-685aa) (CSB-YP3324GMYa4)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human Novel Coronavirus (SARS-CoV-2/ 2019-nCoV)
<b>Tested Applications</b>	ELISA; Recommended dilution: ELISA:1:1000-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	S, S1, S1-RBD, Spike glycoprotein
<b>Immunogen Species</b>	Human Novel Coronavirus (SARS-CoV-2/ 2019-nCoV)
<b>Research Area</b>	Microbiology
<b>Target Names</b>	S
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.