




Recombinant Human Runt-related transcription factor 1 (Runx1)

Product Code	CSB-EP020592HU
Abbreviation	Recombinant Human Runx1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q01196
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MRIPVDASTSRRFTPPSTALSPGKMSEALPLGAPDAGAALAGKLRSGDRSMV EVLADHPGELVRTDSPNFLCSVLPTHWRCNKTLPIAFKVVALGDVPDGLTVTV MAGNDENYSAELRNATAAMKNQVARFNDLRFVGRSGRGKSFTLTITVFTNPP QVATYHRAIKITVDGPREPRRHRQKLLDDQTKPGSLSFSERLSELEQLRRTAMR VSPHHAPTTPNPRASLNHSTAFNPQPQSQMQDTRQIQPSPPWSYDQSYQYL GSIASPSVHPATPISPGRASGMTTSLAELSSRLSTAPDLTAFSDPRQFPALPSIS DPRMHYPGAFTYSPTPVTSGIGIGMSAMGSATRYHTYLPPPYPGSSQAQGGP FQASSPSYHLYYGASAGSYQFSMVGGERSPPRILPPCTNASTGSALLNPSPN QSDVVEAEGSHSNSPTNMAPSARLEEAVWRPY
Research Area	Cancer
Source	E.coli
Target Names	RUNX1
Expression Region	1-453aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	53.8 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.