




ALDH2 Recombinant Monoclonal Antibody

Product Code	CSB-RA799763A0HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05091
Immunogen	A synthesized peptide derived from human ALDH2
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:50-1:200
Relevance	extracellular exosome, mitochondrial matrix, aldehyde dehydrogenase (NAD) activity, aldehyde dehydrogenase [NAD(P)+] activity, electron carrier activity, NAD binding, alcohol metabolic process, carbohydrate metabolic process, ethanol oxidation
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Rabbit IgG in 10mM phosphate buffered saline , pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.
Purification Method	Affinity-chromatography
Isotype	Rabbit IgG
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer; Cardiovascular; Tags & Cell Markers; Metabolism; Signal transduction
Target Names	ALDH2
Clone No.	3F10

Image



Western Blot

Positive WB detected in: Ntera-2 whole cell lysate, Mouse Liver whole cell lysate, Rat Liver cell lysate


All lanes: ALDH2 antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 57, 51 kDa

Observed band size: 57 kDa



IHC image of CSB-RA799763A0HU diluted at 1:100 and staining in paraffin-embedded human liver cancer performed on a Leica Bond™ system. After dewaxing and hydration, antigen retrieval was mediated by high pressure in a citrate buffer (pH 6.0). Section was blocked with 10% normal goat serum 30min at RT. Then primary antibody (1% BSA) was incubated at 4? overnight. The primary is detected by a Goat anti-rabbit IgG polymer labeled by HRP and visualized using 0.05% DAB.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.