




MYC Monoclonal Antibody

Product Code	CSB-MA015270A0m
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01106
Immunogen	Recombinant Human MYC protein (1-439AA)
Raised In	Mouse
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB: 1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG1
Clonality	Monoclonal
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MYC
Clone No.	5H5G5

Image



Western Blot

Positive WB detected in: HeLa whole cell lysate
 All lanes: MYC antibody at 1:500
 Secondary
 Goat polyclonal to mouse IgG at 1/50000 dilution
 Predicted band size: 49, 51 KDa
 Observed band size: 49 KDa
 Exposure time:5min

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.