




Recombinant Pantoea ananas Zeaxanthin glucosyltransferase (crtX)

Product Code	CSB-BP325132ERA
Abbreviation	Recombinant Pantoea ananas crtX protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P21686
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Protein
Immunogen Species	Pantoea ananas (Erwinia uredovora)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MSHFAAIAPPFYSHVRALQNLAQELVARGHRVTFIQQYDIKHLIDSETIGFHSVG TDSHPPGALTRVLHLAAHPLGPSMLKLINEMARTTDMLCRELPQAFNDLAVDG VVDQMEPAGALVAEALGLPFISVACALPLNREPDMP LAVMPFEYGTSDAARE RYAASEKIYDWLMRRHDRVIAEHS HRMGLAPRQKLHQCF SPLAQISQLVPELD FPRKALPACFHAVGPLRETHAPSTSSSR YFTSSEKPRIFASLGT LQGHRYGLF KTIVKACEEIDGQLLLAHCGRLTDSQCEELARSRHTQVVDFADQSAALSQAQL AITHGGMNTVLD AINYRTPL LALPLAFDQPGVASRIVYHGIGKRASRFTTSHALA RQMRSLLTNVDFQQRMAKIQTALRLAGGTMAAAD IIEQVMCTGQPVLSGSGY ATAL
Research Area	Others
Source	Baculovirus
Target Names	crtX
Expression Region	1-431aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	51.1 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.