




Recombinant Arabidopsis thaliana UDP-glycosyltransferase 72E2 (UGT72E2)

Product Code	CSB-BP873379DOA
Abbreviation	Recombinant Mouse-ear cress UGT72E2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LVR1
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MHITKPHAAMFSSPGMGHVIPVIELGKRLSANNGFHVTVFVLETDAASAQSKFL NSTGVDIVKLPSPDIYGLVDPDDHVVTKIGVIMRAAVPALRSKIAAMHQKPTALI VDLFGTDALCLAKEFNMLS YVFIPTNARFLG VSIYYPNLDKDIKEEHTVQRNPLA IPGCEPVR FEDTL DAYLVPDEPVYRDFVRHGLAYPKADGILVNTWEEMEPKSL KSLLNPKLLGRVARVPVYPIGPLCRPIQSSETDHPVLDWLNEQPNESVLYISFG SGGCLSAKQLTELAWGLEQSQR FVWVVRPPVDGSCCSEYVSANGGGTED NTPEYLPEGFVSRTSDRGFV VPSWAPQAEILSHRAVGGFLTHCGWSSTLESV VGGVPMIAWPLFAEQNMNAALLSDELGIAVRLDDPKEDISRWKIEALVRKVMT EKEGEAMRRKVKKLRDSAEMSLSIDGGGLAHESLCRVTKCQRFLERVVDLS RGA
Research Area	Others
Source	Baculovirus
Target Names	UGT72E2
Expression Region	1-481aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	56.9 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.