




Recombinant Mouse Serine/threonine-protein kinase PLK2 (Plk2)

Product Code	CSB-CF018194MO
Abbreviation	Recombinant Mouse Plk2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53351
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MELLRTITYQPAAGTKMCEQALGKACGGDSKKKRPQQPSEDGQPQAQVTPA APHHHHHHSHSGPEISRIIVDPTTGKRYCRGKVLGKGGFAKCYEMTDLTNNKV YAAKIIPHSRVAKPHQREKIDKEIELHRLHHKHVVQFYHYFEDKENIYILLEYCS RRSMAHILKARKVLTEPEVRYYLRLQIVSGLKYLHEQEILHRDLKLGNNFFINEAME LKVGDFGLAARLEPLEHRRRTICGTPNYLSPEVLNKQGHGCESDIWALGCVMY TMLLGRPPFETTNLKETYRCIREARYTMPSSLLAPAKHLIASMLSKNPEDRPSL DDIIRHDFFLQGFTPDLRLSSSCCHTVPDFHLSSPAKNFFFKAAAAALFGGKKDKA RYNDTHNKVSKEDEDIYKLRHDLKKVSITQQPSKHLADEEPQPPPTTVARSGT SAVENKQQIGDAIRMIVRGTLGSCSSSSECLESTMGSVADTVARVLRGCLN MPEADCIPKEQLSTSFQWVTKWVDYSNKYGFYQLSDHTVGVLFNNGAHMS LLPDKKTVHYAELGQCSVFPATDAPEQFISQVTVLKYFSHYMEENLMDGGDL PSVTDIRRPRLYLLQWLKSDKALMMLFNDGTFQVNFYHDHTKIIICNQSEEYLL TYINEDRISTTFRLTLLMSGCSLELKNRMEYALNMLLQRCN
Research Area	Signal Transduction
Source	in vitro E.coli expression system
Target Names	Plk2
Expression Region	1-682aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	85.3 kDa
Protein Length	Full Length



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.