




ACSL6 Antibody

Product Code	CSB-PA892139LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKU0
Immunogen	Recombinant Human Long-chain-fatty-acid--CoA ligase 6 protein (214-484aa)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ACSL6

Image



Western Blot

Positive WB detected in: K562 whole cell lysate
 All lanes: ACSL6 antibody at 1:1000
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 78, 81, 75, 77, 80, 70 kDa
 Observed band size: 78 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.