




# ATAD3A Antibody

<b>Product Code</b>	CSB-PA878897LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9NVI7
<b>Immunogen</b>	Recombinant Human ATPase family aaA domain-containing protein 3A protein (2-240aa)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Product Type</b>	Polyclonal
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	ATAD3A

## Image



### Western Blot

Positive WB detected in: 293T whole cell lysate,

A549 whole cell lysate, K562 whole cell lysate

All lanes: ATAD3A antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 72, 67, 58 kDa

Observed band size: 72 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.