



Avpr1b Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA002469LD01RA |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P48974 |
| Immunogen | Recombinant Rat Vasopressin V1b receptor protein (224-284aa) |
| Raised In | Rabbit |
| Species Reactivity | Rat |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Product Type | Polyclonal |
| Immunogen Species | Rattus norvegicus (Rat) |
| Research Area | Others |
| Target Names | Avpr1b |