



Avpr1b Antibody, FITC conjugated

Product Code	CSB-PA002469LC01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48974
Immunogen	Recombinant Rat Vasopressin V1b receptor protein (224-284aa)
Raised In	Rabbit
Species Reactivity	Rat
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Others
Target Names	Avpr1b
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.