



C1QTNF9 Antibody, Biotin conjugated

Product Code	CSB-PA003654PD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C862
Immunogen	Recombinant Human Complement C1q and tumor necrosis factor-related protein 9A protein (58-72aa)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	C1QTNF9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.