



C22orf39 Antibody, Biotin conjugated

Product Code	CSB-PA761294LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P5X5
Immunogen	Recombinant Human UPF0545 protein C22orf39 protein (1-142aa)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	C22orf39
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.