



C3orf38 Antibody, Biotin conjugated

Product Code	CSB-PA708479LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5JPI3
Immunogen	Recombinant Human Uncharacterized protein C3orf38 protein (80-329aa)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	C3orf38
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.