



CYP2A13 Antibody, FITC conjugated

Product Code	CSB-PA621969LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16696
Immunogen	Recombinant Human Cytochrome P450 2A13 protein (252-354aa)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	CYP2A13
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.