




# ERVFRD-1 Antibody

<b>Product Code</b>	CSB-PA352689LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P60508
<b>Immunogen</b>	Recombinant Human Syncytin-2 protein (33-315aa)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Product Type</b>	Polyclonal
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Others
<b>Target Names</b>	ERVFRD-1

## Image



### Western Blot

Positive WB detected in: HepG2 whole cell lysate, SH-SY5Y whole cell lysate, HT-29 whole cell lysate, 293T whole cell lysate, K562 whole cell lysate, Mouse liver tissue

All lanes: ERVFRD-1 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 60 kDa

Observed band size: 60 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.