



eutC Antibody, Biotin conjugated

Product Code	CSB-PA323953LD01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19636
Immunogen	Recombinant Escherichia coli Ethanolamine ammonia-lyase light chain protein (1-295aa)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Escherichia coli(strain K12)
Research Area	Others
Target Names	eutC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.