



FKBP1B Antibody, HRP conjugated

Product Code	CSB-PA008692LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68106
Immunogen	Recombinant Human Peptidyl-prolyl cis-trans isomerase FKBP1B protein (1-108aa)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	FKBP1B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.