



Fn1 Antibody, Biotin conjugated

Product Code	CSB-PA008759HD01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04937
Immunogen	Recombinant Rat Fibronectin protein (1175-1451aa)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Others
Target Names	Fn1