



Fn1 Antibody, HRP conjugated

Product Code	CSB-PA008759HB01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04937
Immunogen	Recombinant Rat Fibronectin protein (1175-1451aa)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Others
Target Names	Fn1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.