



H6PD Antibody, HRP conjugated

Product Code	CSB-PA010111LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95479
Immunogen	Recombinant Human GDH/6PGL endoplasmic bifunctional protein (422-712aa)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	H6PD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.