




hns Antibody

Product Code	CSB-PA359764XA11ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0ACF8
Immunogen	Recombinant Escherichia coli DNA-binding protein H-NS protein (2-137aa)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Escherichia coli(strain K12)
Research Area	Others
Target Names	hns

Image



Western Blot

Positive WB detected in: Escherichia coli lysate
All lanes: HNS antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 16 kDa
Observed band size: 16 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.