



KBTBD7 Antibody, HRP conjugated

Product Code	CSB-PA823893LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WVZ9
Immunogen	Recombinant Human Kelch repeat and BTB domain-containing protein 7 protein (475-684aa)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	KBTBD7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.