



KCNJ12 Antibody, FITC conjugated

Product Code	CSB-PA619877LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14500
Immunogen	Recombinant Human ATP-sensitive inward rectifier potassium channel 12 protein (182-433aa)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNJ12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.