



KMT5C Antibody, Biotin conjugated

Product Code	CSB-PA769808LD21HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86Y97
Immunogen	Recombinant Human Histone-lysine N-methyltransferase KMT5C protein (235-345aa)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	KMT5C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.