



LCMT1 Antibody, FITC conjugated

Product Code	CSB-PA012808LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UIC8
Immunogen	Recombinant Human Leucine carboxyl methyltransferase 1 protein (157-357aa)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LCMT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.