



MSRB2 Antibody, FITC conjugated

Product Code	CSB-PA896505HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3D2
Immunogen	Recombinant Human Methionine-R-sulfoxide reductase B2, mitochondrial protein (21-182aa)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MSRB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.