



OR5H6 Antibody, FITC conjugated

Product Code	CSB-PA822807PC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGV6
Immunogen	Recombinant Human Olfactory receptor 5H6 protein (277-292aa)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	OR5H6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.