



OR9G1 Antibody, Biotin conjugated

Product Code	CSB-PA818815OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NH87
Immunogen	Recombinant Human Olfactory receptor 9G1 protein (290-304aa)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	OR9G1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.