



OSTN Antibody, FITC conjugated

Product Code	CSB-PA017269LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61366
Immunogen	Recombinant Human Osteocrin protein (28-133aa)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	OSTN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.