



# otsB Antibody, HRP conjugated

Product Code	CSB-PA329093XB01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31678
Immunogen	Recombinant Escherichia coli Trehalose-6-phosphate phosphatase protein (1-266aa)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Escherichia coli(strain K12)
Research Area	Others
Target Names	otsB