



SCNN1D Antibody, FITC conjugated

Product Code	CSB-PA020850LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51172
Immunogen	Recombinant Human Amiloride-sensitive sodium channel subunit delta protein (292-519aa)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SCNN1D
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.