



ZP2 Antibody, Biotin conjugated

Product Code	CSB-PA29809D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05996
Immunogen	Recombinant Human Zona pellucida sperm-binding protein 2 protein (440-716aa)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.