



fbpA Antibody

Product Code	CSB-PA182794ZA01GKB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WQP3
Immunogen	Recombinant Mycobacterium tuberculosis fbpA protein
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Mycobacterium tuberculosis (strain CDC 1551 / Oshkosh)
Lead Time	Made-to-order (12-14 weeks)
Target Names	fbpA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.