



PTP3 Antibody

Product Code	CSB-PA125444ZA01NHN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	R0KX08
Immunogen	Recombinant Nosema bombycis PTP3 protein
Raised In	Rabbit
Species Reactivity	Nosema bombycis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Nosema bombycis (strain CQ1 / CVCC 102059) (Microsporidian parasite)
Lead Time	Made-to-order (12-14 weeks)
Target Names	PTP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.