



Structural polyprotein Antibody

Product Code	CSB-PA182094ZA01ASG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15480
Immunogen	Recombinant Avian infectious bursal disease virus Structural polyprotein protein
Raised In	Rabbit
Species Reactivity	Avian infectious bursal disease virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Avian infectious bursal disease virus (strain Cu-1) (IBDV) (Gumboro disease virus)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Structural polyprotein,partial
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.