



D-lactate dehydrogenase Antibody

Product Code	CSB-PA186974ZA01LPH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51011
Immunogen	Recombinant Leuconostoc mesenteroides subsp. cremoris D-lactate dehydrogenase protein
Raised In	Rabbit
Species Reactivity	Leuconostoc mesenteroides subsp. cremoris
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Leuconostoc mesenteroides subsp. cremoris
Lead Time	Made-to-order (12-14 weeks)
Target Names	Cremoris D-lactate dehydrogenase
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.