



Uncharacterized Antibody

Product Code	CSB-PA187544ZA01CQG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8FNM6
Immunogen	Recombinant Corynebacterium efficiens Uncharacterized protein
Raised In	Rabbit
Species Reactivity	Corynebacterium efficiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Corynebacterium efficiens (strain DSM 44549 / YS-314 / AJ 12310 / JCM 11189 / NBRC 100395)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	Uncharacterized protein