



Capsid Antibody

Product Code	CSB-PA365998ZA01FSB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03612
Immunogen	Recombinant Escherichia phage MS2 Capsid protein
Raised In	Rabbit
Species Reactivity	Escherichia phage MS2
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Escherichia phage MS2 (Bacteriophage MS2)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Capsid protein
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.