



IL17A Antibody

Product Code	CSB-PA011597ZA01MOW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	F6T3G5
Immunogen	Recombinant Macaca mulatta IL17A protein
Raised In	Rabbit
Species Reactivity	Macaca mulatta
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Macaca mulatta (Rhesus macaque)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	IL17A