



APOE Antibody

Product Code	CSB-PA001936ZA01RB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18287
Immunogen	Recombinant Oryctolagus cuniculus APOE protein
Raised In	Rabbit
Species Reactivity	Oryctolagus cuniculus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	APOE