




Barstar (Ribonuclease inhibitor) Antibody

Product Code	CSB-PA320008ZA01BQB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11540
Immunogen	Recombinant Bacillus amyloliquefaciens Barstar (Ribonuclease inhibitor) protein (2-90aa)
Species Reactivity	Bacillus amyloliquefaciens
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bacillus amyloliquefaciens (Bacillus velezensis)
Research Area	Others
Target Names	N/A

Image



Western Blot

Positive WB detected in Recombinant protein
All lanes: Barstar (Ribonuclease inhibitor)
antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 16 kDa
Observed band size: 16 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.