



GPC Antibody

Product Code	CSB-PA362127ZA01LKZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07399
Immunogen	Recombinant Lymphocytic choriomeningitis virus GPC protein
Raised In	Rabbit
Species Reactivity	Lymphocytic choriomeningitis virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Lymphocytic choriomeningitis virus (strain WE) (LCMV)
Lead Time	Made-to-order (12-14 weeks)
Target Names	GPC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.