




Phospholipase A2 Antibody

Product Code	CSB-PA018091ZA01DNK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00630
Immunogen	Recombinant Apis mellifera Phospholipase A2 protein (34-167)
Raised In	Rabbit
Species Reactivity	Apis mellifera
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Apis mellifera (Honeybee)
Target Names	Phospholipase A2

Image



Western Blot

Positive WB detected in: recombinant protein
 All lanes: Phospholipase A2 Antibody at 1:1000
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 20 kDa
 Observed band size: 20 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.