



all1616 Antibody

Product Code	CSB-PA29648ZA01FUZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8YWJ7
Immunogen	Recombinant Nostoc sp. all1616 protein
Raised In	Rabbit
Species Reactivity	Nostoc sp.
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Nostoc sp. (strain PCC 7120 / SAG 25.82 / UTEX 2576)
Lead Time	Made-to-order (12-14 weeks)
Target Names	all1616
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.