




Genome polyprotein Antibody

Product Code	CSB-PA312464ZA01BKX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01499
Immunogen	Recombinant Bovine viral diarrhea virus Genome polyprotein protein (1-168aa)
Species Reactivity	Bovine viral diarrhea virus
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bovine viral diarrhea virus (strain SD-1) (BVDV) (Mucosal disease virus)
Research Area	Others
Target Names	N/A

Image



Western Blot

Positive WB detected in Recombinant protein
All lanes: Genome polyprotein antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 29 kDa
Observed band size: 29 kDa