



HWP1 Antibody

Product Code	CSB-PA343212ZA01CZD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46593
Immunogen	Recombinant Candida albicans HWP1 protein
Raised In	Rabbit
Species Reactivity	Candida albicans
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Candida albicans (strain SC5314 / ATCC MYA-2876) (Yeast)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	HWP1