




A33R Antibody

Product Code	CSB-PA300755ZA01VAA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68616
Immunogen	Recombinant Vaccinia virus (strain Copenhagen) (VACV) A33R protein (57-185aa)
Raised In	Rabbit
Species Reactivity	Vaccinia virus (strain Copenhagen) (VACV)
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Vaccinia virus (strain Copenhagen) (VACV)
Target Names	A33R Antibody

Image



Western Blot

Positive WB detected in: recombinant protein

All lanes:A33R Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 23 kDa

Observed band size: 23 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.