



Hyaluronidase Antibody

Product Code	CSB-PA308534ZA01VBQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P86875
Immunogen	Recombinant Vespa magnifica Hyaluronidase protein
Raised In	Rabbit
Species Reactivity	Vespa magnifica
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Vespa magnifica (Hornet)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Hyaluronidase
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.